




Presentation Part Number: WM 225.163

Medical Combo arm Standard adjustable in height by 300 mm, Mounting on Amico medical device pendants and headwalls. Supported weight balancing range from 6 to 10 Kg





This 300 mm long height adjustable Combo Medical Arm allows to fix, orientate and adjust perfectly a screen, its keyboard support with dimensions of $500 \times 202 \text{ mm}$ delivered with this configuration allows to welcome a standard keyboard as well as a mouse. As an option, it can be equipped with a stand with a sliding right/left mouse tray with dimensions of $480 \times 220 \text{ mm}$ for the keyboard and $226 \times 220 \text{ mm}$ for the mouse, or with a large keyboard stand with dimensions of $580 \times 200 \text{ mm}$. It has a safety stop button for the height adjustment, which allows the arm to remain in the selected position.

Technical specifications:

* Integrated cable management * Standard keyboard support of 500 x 202 mm (See Option on the following pages) * This standard Combo Arm with height adjustment of 300 mm in length simply attaches to the Rail with the Amico adapter and then adjusts and fixes.* The head of this arm conforms to the VESA 75/100mm screen mounting standard * Compatible with equipment weighing up to 6 - 10 Kg * Screen tilt: 10° * Arm rotation: 90° to the right, 90° to the left from the mounting point * Keyboard and screen support rotation: 100° to the right, 100° to the left * Total length of the configuration 857 mm * Colour RAL 5013 cobalt blue and RAL 9016 traffic white * Warranty: 5 years









Technical specifications

Medical Combo arm Standard adjustable in height by 300 mm, Mounting on Amico medical device pendants and headwalls. Supported weight balancing range from 6 to 10 Kg









Monitor adaptation: VESA 75/100 mm Colour: Decorative parts: RAL 5013 cobalt blue Aluminium parts: RAL 9016 traffic white Balance scale of this medical arm: 6 - 10 Kg

Product weight: 4,9 - 6,5 Kg (depending on the dimensions of the keyboard tray)

Height adjustment: Parallel mechanism for a constant viewing angle

Safety: Has a safety stop button for height adjustment, which allows the medical arm to

remain in the chosen position.











Mounting Information - Adapter Mount Amico

Medical Combo arm Standard adjustable in height by 300 mm, Mounting on Amico medical device pendants and headwalls. Supported weight balancing range from 6 to 10 Kg

Adapter Amico:

Mounting on Amico medical device pendants and headwalls







Product highlights:

The height adjustable arms offer full cable integration and easy cleaning for infection control in all critical care areas of the clinical environment. Our 5-year warranty ensures reliable, trouble-free and maintenance-free operation. The swivel function at the point of attachment as well as the tilt and swivel function of the monitor adaptation offer optimal positioning of your device. The height-adjustable arm in its extended version is a two-part arm with a central pivot point between the two arm sections. All bottom mounted monitors have the safety feature of a 20° tilt stop.







Mounting Information - Mounting Options

Medical Combo arm Standard adjustable in height by 300 mm, Mounting on Amico medical device pendants and headwalls. Supported weight balancing range from 6 to 10 Kg

Where do you want to mount the support arm?

At the hospital, there are different ways to fix our brackets: round tube, standard rail, wall rail or on an anaesthesia machine. We will always provide you with the adapter you need.



Adapter for vertical wall slide For wall slide mounting



Adapter for Amico remote controls and wall panels



Adapter for Mindray medical power supply units






Fixed clamp for round tube



For mounting on round tube ø 25, 35, 38 mm







Adapter for Maquet anaesthesia machines Cone Adapter for Model FLOW i



Adapter for Löwenstein anaesthesia machines Medical Adapter for Leon / Leon plus



Universal round tube clamp V-clamp for mounting on round tube ø 23 - 40 mm







Mounting Information - Mounting Options

Medical Combo arm Standard adjustable in height by 300 mm, Mounting on Amico medical device pendants and headwalls. Supported weight balancing range from 6 to 10 Kg

Where do you want to mount the support arm?


At the hospital, there are different ways to fix our brackets: round tube, standard rail, wall rail or on an anaesthesia machine. We will always provide you with the adapter you need.



Adapter for anaesthesia machine (all suppliers) H 80 mm, with existing side rail



Clamp for a vertical standard rail For mounting on standard rail 10 x 25 mm





Considered to be the world's best in advanced hygiene control, our medical mounting systems include cable integration, a 5-year maintenance-free warranty and medical grade compliance. We're simply the best in the world at what we do.







Technical drawing of this Combo configuration

Medical Combo arm Standard adjustable in height by 300 mm, Mounting on Amico medical device pendants and headwalls. Supported weight balancing range from 6 to 10 Kg

