




Presentation Part Number: WM 226.107

Medical Combo 2 x 300 mm swivel arm foldable for monitor and keyboard, (weight supported up to 18 Kg), floor or ceiling column mount







Technical specifications:

- * Integrated cable management
- * Standard 500 x 202 mm keyboard support foldable 90°
- * Optional floor or ceiling columns are available in various heights (1500 2000 mm). These columns can be easily cut to size. All cables are also routed through the column. (See Option on the following pages)
- * This 300mm long Combo Arm simply slides into the column rail and then adjusts to the required height and is fixed. It does not have a height adjustment but does have the ability to rotate.
- * The head of this arm complies with the VESA 75/100mm screen mounting standard.
- * Compatible with equipment weighing up to 18kg.
- * Screen tilt: 10° * Arm rotation: 105° right, 105° left. 180° central rotation
- * Keyboard support rotation: 110° right, 110° left
- * Total length of the configuration 986 mm * Colour RAL 5013 cobalt blue and RAL 9016 traffic white
- * Warranty: 5 years










Technical specifications

Medical Combo 2 x 300 mm swivel arm foldable for monitor and keyboard, (weight supported up to 18 Kg), floor or ceiling column mount









Monitor adaptation: VESA 75/100 mm

Colour: Les pièces décoratives: RAL 5013 cobalt blue

Aluminium parts: RAL 9016 traffic white **Max. load capacity of the support arm:** 18 kg **Monitor adaptation Max. load capacity:** 18 Kg

Product weight: 4,7 to 6,3 Kg - depending on the size of the keyboard tray

Product highlights:

Non-height adjustable arm. The extended swivel arms and swivel arms provide full cable integration and easy cleaning for infection control in all ICU areas of the clinical environment. Our 5-year warranty guarantees reliable, trouble-free and maintenance-free operation. The swivel function at the attachment point as well as the tilt and swivel function of the monitor adapter provide optimal positioning of your device. The extended swing arm is a two-part arm with a central pivot point between the two sections of the arm. All bottom-mounted monitors have the safety function of a 20° tilt stop.









Mounting Information - Column Slide Rail

Medical Combo 2 x 300 mm swivel arm foldable for monitor and keyboard, (weight supported up to 18 Kg), floor or ceiling column mount

Floor or ceiling column:


To install this medical arm, you need to use one of our medical columns. These are available in standard heights of 1500 and 2000 mm. Custom heights are available on request. Slide covers and end caps are supplied with every order.

















Mounting Information - Mounting Options

Medical Combo 2 x 300 mm swivel arm foldable for monitor and keyboard, (weight supported up to 18 Kg), floor or ceiling column mount

Where do you want to mount the support arm?

At the hospital, there are different ways to fix our brackets: round tube, standard rail, wall rail or on an anaesthesia machine. We will always provide you with the adapter you need.



Adapter for vertical wall slide
For wall slide mounting



Adapter for Amico remote controls and wall panels



Adapter for Mindray medical power supply units



Double clamp for vertical standard rail

For mounting on two parallel vertical standard rails 10 x 25 mm



Fixed clamp for a standard rail

For mounting on a horizontal standard rail 10 x 25, 8 x 35, 10 x 30, 10 x 50 mm



Fixed clamp for round tube

For mounting on round tube ø 25, 35, 38 mm



Adapter for Maquet anaesthesia machines Cone Adapter for Model FLOW i



Adapter for Löwenstein anaesthesia machines Medical Adapter for Leon /

Medical Adapter for Leon Leon plus



Universal round tube clamp

V-clamp for mounting on round tube ø 23 - 40 mm







Mounting Information - Mounting Options

Medical Combo 2 x 300 mm swivel arm foldable for monitor and keyboard, (weight supported up to 18 Kg), floor or ceiling column mount

Where do you want to mount the support arm?


At the hospital, there are different ways to fix our brackets: round tube, standard rail, wall rail or on an anaesthesia machine. We will always provide you with the adapter you need.



Adapter for anaesthesia machine (all suppliers) H 80 mm, with existing side rail



Clamp for a vertical standard rail For mounting on standard rail 10 x 25 mm





Considered to be the world's best in advanced hygiene control, our medical mounting systems include cable integration, a 5-year maintenance-free warranty and medical grade compliance. We're simply the best in the world at what we do.







Technical drawing of this Combo configuration

Medical Combo 2 x 300 mm swivel arm foldable for monitor and keyboard, (weight supported up to 18 Kg), floor or ceiling column mount

