



Presentation Part Number: WM 226.111

Medical Combo arm Foldable and adjustable in height by 300 mm, floor or ceiling column mounting. Supported weight of balance scale from 13,5 to 16,5 Kg





This 300 mm long height adjustable and foldable Combo Medical Arm allows to fix, orientate and adjust perfectly a screen, its keyboard support with dimensions of 500 x 202 mm delivered with this configuration allows to welcome a standard keyboard as well as a mouse. As an option, it can be equipped with a stand with a sliding right/left mouse tray with dimensions of 480 x 220 mm for the keyboard and 226 x 220 mm for the mouse, or with a large keyboard stand with dimensions of 580 x 200 mm. It has a safety stop button for the height adjustment, which allows the arm to remain in the selected **position**.

Technical specifications:

* Integrated cable management * Standard keyboard support 500 x 202 mm foldable 90° * Optional floor or ceiling columns are available in various heights (1500 - 2000 mm). These columns can be easily cut to size. All cables are also routed through the column. (See Option on the following pages) * This folding Combo Arm with height adjustment is 300 mm long and simply slides into the column rail and then adjusts to the required height and is fixed * The head of this arm complies with the screen mounting standard: VESA 75/100 mm * Compatible with equipment weighing up to 13,5 - 16,5 kg. * Tilting of the screen: 10° * Rotation of the arm: 90° to the right, 90° to the left from the fixing point * Rotation of the keyboard and screen support: 100° to the right, 100° to the left * Total length of the configuration 977 mm * Colour RAL 5013 cobalt blue and RAL 9016 traffic white * Warranty: 5 years









Technical specifications

Medical Combo arm Foldable and adjustable in height by 300 mm, floor or ceiling column mounting. Supported weight of balance scale from 13,5 to 16,5 Kg



Cette colonne possède deux rails de fixation pour les différents éléments qui y seront fixés. (1 rail à l'avant de la colonne, 1 rail à l'arrière) This column has two fixing rails for the different elements that will be fixed to it (1 rail at the front of the column, 1 rail at the back) Deze zuil heeft twee bevestigingsrails voor de verschillende elementen die eraan bevestigd zullen worden (1 rail aan de voorkant van de zuil, 1 rail aan de achterkant)





Monitor adaptation: VESA 75/100 mm

Colour: Decorative parts: RAL 5013 cobalt blue

Aluminium parts: RAL 9016 traffic white

Balance scale of this medical arm: 13,5 - 16,5 Kg

Product weight: 4,9 - 6,5 Kg (depending on the dimensions of the keyboard tray)

Height adjustment: Parallel mechanism for a constant viewing angle Safety: Has a safety stop button for height adjustment, which allows the medical arm to remain in the chosen position.













Mounting Information - Column Slide Rail

Medical Combo arm Foldable and adjustable in height by 300 mm, floor or ceiling column mounting. Supported weight of balance scale from 13,5 to 16,5 Kg

Floor or ceiling column:


To install this medical arm, you need to use one of our medical columns. These are available in standard heights of 1500 and 2000 mm. Custom heights are available on request. Slide covers and end caps are supplied with every order.

















Mounting Information - Mounting Options

Medical Combo arm Foldable and adjustable in height by 300 mm, floor or ceiling column mounting. Supported weight of balance scale from 13,5 to 16,5 Kg

Where do you want to mount the support arm?

At the hospital, there are different ways to fix our brackets: round tube, standard rail, wall rail or on an anaesthesia machine. We will always provide you with the adapter you need.



Adapter for vertical wall slide For wall slide mounting



Adapter for Amico remote controls and wall panels



Adapter for Mindray medical power supply units









Double clamp for vertical standard rail For mounting on two parallel vertical standard

rails 10 x 25 mm



For mounting on a horizontal standard rail 10 x 25, 8 x 35, 10 x 30, 10 x 50 mm



tube ø 25, 35, 38 mm



Universal round tube clamp V-clamp for mounting on

round tube ø 23 - 40 mm



Adapter for Maquet anaesthesia machines Cone Adapter for Model FLOW i



Adapter for Löwenstein anaesthesia machines Medical Adapter for Leon / Leon plus





Mounting Information - Mounting Options

Medical Combo arm Foldable and adjustable in height by 300 mm, floor or ceiling column mounting. Supported weight of balance scale from 13,5 to 16,5 Kg

Where do you want to mount the support arm?

At the hospital, there are different ways to fix our brackets: round tube, standard rail, wall rail or on an anaesthesia machine. We will always provide you with the adapter you need.



Adapter for anaesthesia machine (all suppliers) H 80 mm, with existing side rail



Clamp for a vertical standard rail For mounting on standard rail 10 x 25 mm



Product Highlights:

Height adjustable arms offer full cable integration and easy cleaning for infection control in all critical care areas of the clinical environment. Our 5-year warranty ensures reliable, trouble-free and maintenance-free operation. The swivel function at the point of attachment as well as the tilt and swivel function of the monitor adaptation offer optimal positioning of your device. The height-adjustable arm in its extended version is a twopart arm with a central pivot point between the two arm sections. All bottom mounted monitors have the safety feature of a 20° tilt stop.



Considered to be the world's best in advanced hygiene control, our medical mounting systems include cable integration, a 5-year maintenance-free warranty and medical grade compliance. We're simply the best in the world at what we do.







Technical drawing of this Combo configuration

Medical Combo arm Foldable and adjustable in height by 300 mm, floor or ceiling column mounting. Supported weight of balance scale from 13,5 to 16,5 Kg

